



AESCULAP® POWER SYSTEMS

CATALOG 2020

ELAN 4 and Acculan 4

CONTENT

ELAN 4

- ELAN 4 electro
- ELAN 4 air
- ELAN 4 Accessories

Acculan 4

- Acculan 4 Et Acculan 4 Small Drill
- Acculan 4 Electro
- Acculan 4 Accessories

PROCESSING

- Sterilit® Power Systems
- Processing

OTHERS

- Tool Attachments Overview
 - Signs Et Symbols
 - Brochures
 - Index
-

AESCULAP® ELAN 4

PRECISION. PERFORMANCE. PLEASURE.

ELAN 4 is the next generation of pneumatic and electric power systems for neuro and spine departments. With a new technical concept ELAN 4 addresses important requirements of the surgeon and all parties involved in handling power systems in the hospital. ELAN 4 electro and ELAN 4 air are two systems with one handling philosophy.

- Convenient and reliable grip
- Powerful, precise and reliable work
- Easy to set up and simple to use



AESCULAP® ELAN 4

ELAN 4 electro

ELECTRIC HIGHSPEED SYSTEM



ELAN 4 electro

BASIC COMPONENTS

GA800

ELAN 4 electro CONTROL UNIT

with irrigation pump; touch screen;
power cord not included
pump flow 0–65 ml/min

100–240 V AC; 50/60 Hz;
330 x 330 x 200 mm; 9,500 g

consisting parts:
GD412804 bottle holder



GA808

ELAN 4 electro FOOT CONTROL

protection type IPX8; anaesthetic
proof class AP

5.0 m;
245 x 144 x 76 mm; 2,240 g



GA810

ELAN 4 electro WIRELESS FOOT CONTROL

for use with ELAN 4 electro applied
parts

271 x 266 x 85 mm; 3,800 g

Requires control unit with software
version V3.00.



GA806

GB073R

ELAN 4 electro MOTOR CABLE

for use with ELAN 4 electro applied
parts

4.0 m; 270 g

consisting parts:
GB600860 Sterilit® spray adaptor
ELAN 4 electro



AESCULAP® ELAN 4

ELAN 4 electro

BASIC COMPONENTS, PERFORATOR DRIVER & LOWSPEED MOTOR

GA804

GB073R

ELAN 4 electro MOTOR CABLE WITH HAND SWITCH

for use with ELAN 4 electro applied parts

4.0 m; 370 g

consisting parts:
GB600860 Sterilit® spray adaptor
ELAN 4 electro



GA805

GB073R

ELAN 4 electro MOTOR CABLE WITH HAND CONTROL

for use with ELAN 4 electro applied parts

4.0 m; 370 g

consisting parts:
GB600860 Sterilit® spray adaptor
ELAN 4 electro



GA822

GB076R

ELAN 4 electro PERFORATOR DRIVER

for use with cranial perforators with Hudson shank

180 W; 0 - 1,200 rpm; 250 Ncm;
166 x 35 mm; 410 g



GA824

GB077R

ELAN 4 electro LOWSPEED MOTOR

with intra coupling

180 W; 0 - 40,000* rpm;
10 Ncm; 123 x 23 mm; 150 g

* Requires control unit with software version V3.00. With lower versions the motor runs max. 20,000 rpm.



ELAN 4 electro




CRANIOTOME & ATTACHMENTS

GA849

GB085R

ELAN 4 electro CRANIOTOME AND MULTIFUNCTION HANDPIECE

2-ring coding; for use with dura guards and holding sleeves GB941R – GB947R; spray nozzle GA259SU (not included)

140 W;  ; 0 – 80,000 rpm; 2.2 Ncm; 156 x 17 mm;  121 g



GB941R

GB719R

ELAN 4 DURA GUARD PEDIATRIC

fixed; 1 dot coding; for craniotome and multifunction handpieces GA749, GA849



GB942R

GB719R

ELAN 4 DURA GUARD STANDARD

fixed; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849



GB943R

GB719R

ELAN 4 DURA GUARD LONG

fixed; 3 dot coding; for craniotome and multifunction handpieces GA749, GA849




GB945R

GB719R

ELAN 4 HOLDING SLEEVE

for craniotome and multifunction handpieces GA749, GA849



GB947R

GB719R

ELAN 4 DURA GUARD STANDARD

turnable; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849



AESCULAP® ELAN 4




ELAN 4 electro

1-RING HIGHSPEED HANDPIECES

GA861

GB083R

ELAN 4 electro 1-RING HANDPIECE L4




1-ring coding;
spray nozzle GA261SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
152 x 16 mm;  80 g



GA862

GB083R

ELAN 4 electro 1-RING HANDPIECE L7




1-ring coding;
spray nozzle GA262SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
177 x 16 mm;  83 g



GA863

GB083R

ELAN 4 electro 1-RING HANDPIECE L10




1-ring coding;
spray nozzle GA263SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
206 x 16 mm;  87 g



GA864

GB083R

ELAN 4 electro 1-RING HANDPIECE L13

1-ring coding;
spray nozzle GA264SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
237 x 16 mm;  91 g






ELAN 4 electro

2-RING HIGHSPEED HANDPIECES

GA865

GB083R

ELAN 4 electro 2-RING HANDPIECE L7




2-ring coding;
spray nozzle GA265SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
176 x 16 mm;  89 g




GA866

GB083R

ELAN 4 electro 2-RING HANDPIECE L10




2-ring coding;
spray nozzle GA266SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
206 x 16 mm;  95 g




GA867

GB083R

ELAN 4 electro 2-RING HANDPIECE L13




2-ring coding;
spray nozzle GA267SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
236 x 16 mm;  100 g



GA868

GB083R

ELAN 4 electro 2-RING HANDPIECE L17




2-ring coding;
spray nozzle GA268SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
276 x 16 mm;  108 g



GA869

GB083R

ELAN 4 electro 2-RING HANDPIECE L22

2-ring coding;
spray nozzle GA269SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
326 x 16 mm;  118 g



AESCULAP® ELAN 4

ELAN 4 electro

MIS HANDPIECE & HANDPIECE SHAFTS



GA860**GB084R****ELAN 4 electro MIS HANDPIECE**for use with handpiece shafts
GB920R - GB927R140 W;  ; 0 - 80,000 rpm;
2.2 Ncm; 118 x 17 mm;  100 g**GB920R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 STRAIGHT**

spray nozzle GB796SU (not included)

**GB921R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 CURVED**

spray nozzle GB796SU (not included)

**GB922R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 HEAVY CURVED**

spray nozzle GB796SU (not included)

**GB925R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 STRAIGHT**

spray nozzle GB797SU (not included)

**GB926R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 CURVED**

spray nozzle GB797SU (not included)

**GB927R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 HEAVY CURVED**

spray nozzle GB797SU (not included)



AESCULAP® ELAN 4


ELAN 4 electro

SAWS



GA836**GB086R****ELAN 4 electro MICRO SAGITTAL SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 161 x 23 mm;
 \angle 7°; 190 g**GA837****GB087R****ELAN 4 electro MICRO RECIPROCATING SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 184 x 23 mm;
 \updownarrow 1.3 mm; 198 g**GA831****GB088R****ELAN 4 electro SAGITTAL SAW**


spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 129 x 23 mm;
 \angle 4.5°; 220 g**GA832****GB087R****ELAN 4 electro RECIPROCATING SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 194 x 23 mm;
 \updownarrow 2.5 mm; 238 g**GA833****GB088R****ELAN 4 electro TRANSVERSAL SAW**

spray nozzle GA258SU (not included)

180 W; $\swarrow \searrow$; 0 - 20,000 rpm; 10 Ncm;
 130 x 23 mm; \angle 10°; 185 gto be used with saw blade shafts
 GP580R, GP582R, GP586R and GP590R

AESCULAP® ELAN 4

ELAN 4 electro

DRILL



GA844**GB482R****ELAN 4 electro DRILL**

with oscillating and tapping mode

200 W;   3.2 mm;130 x 27 x 122 mm;  800 g

consisting parts:

GA344244 K-wire protection guide



GA344211 rinsing adaptor

GB600880 Sterilit® spray adaptor

Requires control unit with software version V3.00.

**GB870R****GB716R****INTRAMEDULLARY REAMER ATTACHMENT**

with Jacobs chuck 0.5 - 7.4 mm, not for use with acetabulum reamers



0 - 320 rpm; 6.5 Nm;  3.2 mm;105 x 29 mm;  230 g

consisting parts:



key GA031R

**GB871R****GB716R****INTRAMEDULLARY REAMER ATTACHMENT**

with large AO chuck, not for use with acetabulum reamers

0 - 320 rpm; 6.5 Nm;  3.2 mm;88 x 29 mm;  170 g**GB872R****GB716R****INTRAMEDULLARY REAMER ATTACHMENT**

with Hudson/Zimmer chuck, not for use with acetabulum reamers

0 - 320 rpm; 6.5 Nm;  3.2 mm;86 x 25 mm;  150 g

AESCULAP® ELAN 4

ELAN 4 electro ATTACHMENTS

GB881R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm

0 – 1,250 rpm; 2 Nm;  3.2 mm;
89 x 29 mm;  163 g

consisting parts:
key GA031R




GB882R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 4 mm

0 – 1,250 rpm; 2 Nm;  3.2 mm;
76 x 22 mm;  100 g

consisting parts:
key GA062R





GB884R

GB716R

DRILL ATTACHMENT

with small AO chuck

0 – 1,250 rpm; 2 Nm;  2.8 mm;
60 x 22 mm;  66 g





GB886R

GB716R

DRILL ATTACHMENT

with AESCULAP® hexagonal chuck

0 – 1,250 rpm; 2 Nm;  3.2 mm;
57 x 24 mm;  67 g





GB887R

GB716R

DRILL ATTACHMENT

with Dental shank

0 – 1,250 rpm; 2 Nm;  1.8 mm;
56 x 22 mm;  50 g



GB878R**GB716R****ADAPTOR**

for Synthes radiolucent drive 511.300

0 - 1,250 rpm; 2 Nm;

111 x 27 mm; 181 g

**GB894R****GB716R****K-WIRE ATTACHMENT**

for K-wires from 0.6 - 3.2 mm

0 - 1,250 rpm; 3.2 mm;

92 x 27 x 124 mm; 179 g

**GB891R****GB716R****SAGITTAL SAW ATTACHMENT**

for "L-sagittal" sawblades with a maximum usable length of 50 mm

0 - 15,750 rpm; 74 x 33.5 mm;

± 4.5°; 179 g

**GB892R****GB716R****RECIPROCATING SAW ATTACHMENT**

for reciprocating sawblades

0 - 15,750 rpm; 135 x 24.5 mm;

2.5 mm; 166 g



AESCULAP® ELAN 4

ELAN 4 electro



ATTACHMENTS

GB896R

GB716R

SCREW ATTACHMENT

with small AO chuck

0 - 165 rpm; 12 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  235 g





GB897R

GB716R

SCREW ATTACHMENT

with AESCULAP® hexagonal chuck

0 - 165 rpm; 12 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  231 g



AESCULAP® ELAN 4

ELAN 4 air

PNEUMATIC HIGHSPEED SYSTEM



ELAN 4 air WALL HOSE & ADAPTOR

GA702R

ELAN 4 air WALL HOSE

for use with ELAN 4 air wall adaptor
GA710R - GA715R

 3.0 m

GA703R

ELAN 4 air WALL HOSE

for use with ELAN 4 air wall adaptor
GA710R - GA715R

 5.0 m



GA710R

ELAN 4 air WALL ADAPTOR

AESCULAP® Draeger

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



GA712R

ELAN 4 air WALL ADAPTOR

Schrader with diffuser

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



GA713R

ELAN 4 air WALL ADAPTOR

DIN with diffuser

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



GA715R

ELAN 4 air WALL ADAPTOR

Synthes

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



AESCULAP® ELAN 4

ELAN 4 air

FOOT CONTROL, MOTOR HOSE & PERFORATOR DRIVER

GA708

ELAN 4 air FOOT CONTROL

for use with ELAN 4 air wall adaptor GA710R – GA715R; ELAN 4 air wall hose GA702R, GA703R; ELAN 4 air motor hose GA705R – GA707R



GA705R

GB053R

ELAN 4 air MOTOR HOSE

for use with ELAN 4 air applied parts and ELAN 4 air foot control GA708

3.0 m

consisting parts:
GB600850 Sterilit® spray adaptor
ELAN 4 air



GA706R

GB053R

ELAN 4 air MOTOR HOSE

with hand control; for use with ELAN 4 air applied parts and ELAN 4 wall adaptor GA710R – GA715R

3.0 m

consisting parts:
GB600850 Sterilit® spray adaptor
ELAN 4 air



GA707R


GB053R

ELAN 4 air MOTOR HOSE

with hand control; for use with ELAN 4 air applied parts and ELAN 4 wall adaptor GA710R – GA715R

5.0 m

consisting parts:
GB600850 Sterilit® spray adaptor
ELAN 4 air



GA722

GB056R

ELAN 4 air PERFORATOR DRIVER

for use with cranial perforators with Hudson shank

100 W; 0 – 1,200 rpm; 400 Ncm;
155 x 35 mm; 427 g



ELAN 4 air

CRANIOTOME & ATTACHMENTS

GA749

GB065R

ELAN 4 air CRANIOTOME AND MULTIFUNCTION HANDPIECE

2-ring coding; for use with dura guards and holding sleeves GB941R – GB947R

80 W; 0 – 80,000 rpm; 3.5 Ncm;
128 x 18 mm; 108 g




GB941R

GB719R

ELAN 4 DURA GUARD PEDIATRIC

fixed; 1 dot coding;
for craniotome and multifunction
handpieces GA749, GA849




GB942R

GB719R

ELAN 4 DURA GUARD STANDARD

fixed; 2 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



GB943R

GB719R

ELAN 4 DURA GUARD LONG

fixed; 3 dot coding;
for craniotome and multifunction
handpieces GA749, GA849




GB945R

GB719R

ELAN 4 HOLDING SLEEVE

for craniotome and multifunction
handpieces GA749, GA849



GB947R

GB719R

ELAN 4 DURA GUARD STANDARD

turnable; 2 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



AESCULAP® ELAN 4

ELAN 4 air



1-RING HIGHSPEED HANDPIECES

GA761

GB063R

ELAN 4 air 1-RING HANDPIECE L4

1-ring coding

80 W; ; 0–80,000 rpm; 3.5 Ncm;
125 x 18 mm;  68 g





GA762

GB063R

ELAN 4 air 1-RING HANDPIECE L7

1-ring coding

80 W; ; 0–80,000 rpm; 3.5 Ncm;
150 x 18 mm;  71 g





GA763

GB063R

ELAN 4 air 1-RING HANDPIECE L10

1-ring coding

80 W; ; 0–80,000 rpm; 3.5 Ncm;
180 x 18 mm;  75 g





GA764

GB063R

ELAN 4 air 1-RING HANDPIECE L13

1-ring coding

80 W; ; 0–80,000 rpm; 3.5 Ncm;
210 x 18 mm;  79 g



ELAN 4 air



2-RING HIGHSPEED HANDPIECES


GA765

GB063R

ELAN 4 air 2-RING HANDPIECE L7

2-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
150 x 18 mm;  76 g






GA766

GB063R

ELAN 4 air 2-RING HANDPIECE L10

2-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
180 x 18 mm;  82 g





GA767

GB063R

ELAN 4 air 2-RING HANDPIECE L13

2-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
210 x 18 mm;  88 g





GA768

GB063R

ELAN 4 air 2-RING HANDPIECE L17

2-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
250 x 18 mm;  95 g





GA769

GB063R

ELAN 4 air 2-RING HANDPIECE L22

2-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
300 x 18 mm;  105 g



AESCULAP® ELAN 4

ELAN 4 air

MIS HANDPIECE & HANDPIECE SHAFTS



GA760**GB064R****ELAN 4 air MIS HANDPIECE**for use with handpiece shafts
GB920R - GB927R80 W; 0 - 80,000 rpm;
3.5 Ncm; 89 x 18 mm; 85 g**GB920R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 STRAIGHT****GB921R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 CURVED****GB922R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 HEAVY CURVED****GB925R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 STRAIGHT****GB926R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 CURVED****GB927R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 HEAVY CURVED**

AESCULAP® ELAN 4

ELAN 4 air

SAWS & LOWSPEED MOTOR



GA736**GB066R****ELAN 4 air MICRO SAGITTAL SAW**

60 W; 0 – 18,000 rpm; 150 x 23 mm;
 \angle 7°; $\bar{\bar{m}}$ 214 g

**GA737****GB067R****ELAN 4 air MICRO RECIPROCATING SAW**

75 W; 0 – 21,000 rpm; 167 x 23 mm;
 $\bar{\bar{m}}$ 1.3 mm; $\bar{\bar{m}}$ 215 g

**GA731****GB068R****ELAN 4 air SAGITTAL SAW**

75 W; 0 – 21,000 rpm; 116 x 23 mm;
 \angle 4.5°; $\bar{\bar{m}}$ 240 g

**GA732****GB067R****ELAN 4 air RECIPROCATING SAW**

75 W; 0 – 21,000 rpm; 118 x 23 mm;
 $\bar{\bar{m}}$ 2.5 mm; $\bar{\bar{m}}$ 250 g

**GA733****GB068R****ELAN 4 air TRANSVERSAL SAW**

75 W; \searrow ; 0 – 21,000 rpm; 7.6 Ncm;
 115 x 23 mm; \angle 10°; $\bar{\bar{m}}$ 200 g

to be used with saw blade shafts
 GP580R, GP582R, GP586R and GP590R

**GA724****GB057R****ELAN 4 air LOWSPEED MOTOR**

with intra coupling

80 W; \searrow ; 0 – 21,000 rpm; 7.6 Ncm;
 93 x 23 mm; $\bar{\bar{m}}$ 130 g



AESCLAP[®] ELAN 4

ELAN 4 ACCESSORIES

FOR ELAN 4 air & ELAN 4 electro

GB861P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE

for use with 1-ring handpieces



GB863P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE

for use with 2-ring handpieces



GB862P

GB727R

ELAN 4 MIS SLEEVE

for use with 1-ring handpieces



GB866R

ELAN 4 LATERAL BURR GUARD

for use with 1-ring handpieces and
burrs up to 4.5 mm diameter



TE780

POWER CORD

Europe, 1.5 m

TE730

POWER CORD

Europe, 5.0 m

TE734

POWER CORD

UK, 5.0 m

TE735


POWER CORD

USA, Canada, Japan, 3.5 m



GA395SUsterile  

ELAN 4 electro DISPOSABLE TUBE SET

**ELAN 4 SPAY NOZZELS**sterile **GA259SU**

for craniotome attachment

GA261SU

for 1-ring handpiece L4

GA262SU

for 1-ring handpiece L7

GA263SU

for 1-ring handpiece L10

GA264SU

for 1-ring handpiece L13

GA265SU

for 2-ring handpiece L7

GA266SU

for 2-ring handpiece L10

GA267SU

for 2-ring handpiece L13

GA268SU

for 2-ring handpiece L17

GA269SU

for 2-ring handpiece L22

GB796SU

for MIS handpiece shafts L10

GB797SU

for MIS handpiece shafts L13

**GA258SU**

for saws



AESCULAP® Acculan 4

LOOK. FEEL. EXPERIENCE.

Acculan 4 the battery-operated power system for orthopaedics, traumatology, cardiothoracic surgery and dermatology

- Powerful
- Longevity – Exclusion of battery and electronic control from sterilization process
- Easy handling – Titanium housing allows mechanical alkaline reprocessing



AESCULAP® Acculan 4

DRILL & ATTACHMENTS

GA330

GB495R

Acculan 4 DRILL AND REAMER

250 W;   4.0 mm;
207 x 121 x 55 mm;  1,250 g

consisting parts:

GA678 sterile funnel

GA672207 K-wire protection sleeve

GA672211 rinsing adaptor

GB600840 Sterilit® spray adaptor

Acculan

Lid GA675 and NiMH battery GA676 to be ordered separately.

Measurements and weight are including battery and lid.





GB620R

GB496R

Acculan 4 DRILL ATTACHMENT

with Jacobs chuck keyless 0.5 – 7.4 mm

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
133 x 35 mm;  397 g





GB621R

GB496R

Acculan 4 DRILL ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
118 x 35 mm;  300 g

consisting parts:

key GA031R





GB623R

GB496R

Acculan 4 DRILL ATTACHMENT

with small AO chuck

0 – 1,250 rpm; 4.0 Nm;  2.8 mm;
90 x 35 mm;  210 g





ATTACHMENTS

GB639R

GB496R

Acculan 4 DRILL ATTACHMENT

with Jacobs chuck 0.5 – 4.0 mm

0 – 1,250 rpm; 4.0 Nm;  4.0 mm;
105 x 35 mm;  243 g

consisting parts:
key GA062R



GB645R

GB496R

Acculan DRILL ATTACHMENT

for Synthes radiolucent drive 511.300

0 – 900 rpm; 110 x 35 mm;  225 g





GB663R

GB496R

Acculan 4 DRILL ATTACHMENT

with trinkle chuck

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
93 x 35 mm;  217 g





GB664R

GB496R

Acculan 4 DRILL ATTACHMENT

with AESCULAP® hexagonal chuck

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
86 x 35 mm;  225 g





GB665R

GB496R

Acculan 4 DRILL ATTACHMENT

with Hudson/Zimmer chuck

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
100 x 35 mm;  230 g



AESCULAP® Acculan 4



ATTACHMENTS

GB657R

GB496R

Acculan 4 INTRAMEDULLAR ATTACHMENT

with large AO chuck

0–370 rpm; 12.0 Nm;  4.0 mm;
108 x 35 mm;  300 g





GB667R


GB496R

Acculan 4 REAMER ATTACHMENT

with Jacobs chuck 0.5–7.4 mm chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
133 x 35 mm;  280 g

consisting parts:
key GA031R





GB668R

GB496R

Acculan 4 REAMER ATTACHMENT

with large AO chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
108 x 35 mm;  300 g





GB669R

GB496R

Acculan 4 REAMER ATTACHMENT

with Hudson/Zimmer chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
118 x 35 mm;  300 g





GB670R

GB496R

Acculan 4 REAMER ATTACHMENT

with Harris chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
106 x 35 mm;  290 g



GB660R**GB496R****Acculan SAGITTAL SAW ATTACHMENT**

for Acculan rapid action saw blades
with a maximum useable length of
50 mm

«» 0 - 17,000 rpm; \angle 4°47';
89 x 36 x 44 mm; $\bar{\bar{m}}$ 270 g

**GB641R****GB496R****Acculan K-WIRE ATTACHMENT**

for K-wires from 0.6 - 4.0 mm

0 - 1,250 rpm; 93 x 36 x 137 mm;
 $\bar{\bar{m}}$ 285 g



AESCULAP® Acculan 4

OSCILLATING SAWS

GA331

GB495R

Acculan 4 OSCILLATING SAW

250 W; ⇄ 0 – 13,000 rpm; ∠ 4.1°;
211 x 177 x 55 mm; 1,600 g

consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GA336

GB498R

Acculan 4 OSCILLATING SAW STRAIGHT

250 W; ⇄ 0 – 13,000 rpm; ∠ 4.1°;
290 x 55 x 84 mm; 1,480 g

consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GA337

GB498R

Acculan 4 OSCILLATING STERNUM SAW

250 W; ⇄ 0 – 13,000 rpm; ∠ 4.1°;
290 x 55 x 84 mm; 1,480 g

consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



RECIPROCATING SAW

GA334

GB495R

Acculan 4 RECIPROCATING SAW

250 W; ⇄ 0 – 15,000 rpm; ↓ 3.2 mm;
207 x 180 x 55 mm; ⚖ 1,540 g

consisting parts:

GA678 sterile funnel

GB440215 cleaning plate

GA674890 cleaning adaptor

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GB436R

GB686R

PROTECTION GUIDE

for use with GA334 and saw blade
GC761R



AESCULAP® Acculan 4

DERMATOME

GA340

GB498R

Acculan 4 DERMATOME

width keyless adjustable;
thickness adjustable (1 / 10 mm steps);
cutting thickness: 0.2 - 1.2 mm;
cutting width: max. 78 mm;
suitable blades: GB228R
(to be ordered separately)

250 W; ⇄ 0 - 6,500 rpm; ⬆ 3.1 mm;
277 x 113 x 83 mm; ⚖ 1,650 g

consisting parts:
GA678 sterile funnel

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.



GA341

GB498R

Acculan 4 DERMATOME 0.1 MM

width keyless adjustable;
thickness adjustable (1 / 10 mm steps);
cutting thickness: 0.1 - 1.1 mm;
cutting width: max. 78 mm;
suitable blades: GB228R
(to be ordered separately)

250 W; ⇄ 0 - 6,500 rpm; ⬆ 3.1 mm;
277 x 113 x 83 mm; ⚖ 1,650 g

consisting parts:
GA678 sterile funnel

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.



BA720R

GB688R

SKIN GRAFT MESHER

for use in combination with skin graft
carriers BA721 (1:1.5) / BA722 (1:3) /
BA723 (1:6)

150 x 195 x 105 mm; ⚖ 4,600 g

consisting parts:
BA725R cutting roller
BA726R ratchet wrench




AESCULAP® Acculan 4

LOOK. FEEL. EXPERIENCE.

Acculan 4 small drill – the smallest member of the Acculan 4 family can be used for small fragment Osteosynthesis as well as for reamed intramedullary nailing.

Acculan 4



AESCULAP® Acculan 4

SMALL DRILL

DRILL & ATTACHMENTS

GA344

GB482R

Acculan 4 SMALL DRILL

with oscillating and tapping mode;
not for use with acetabulum reamers

85 W;  3.2 mm;
177 x 124 x 55 mm;  840 g

consisting parts:

GA678 sterile funnel

GA344244 K-wire protection guide

GA344211 rinsing adaptor

GB600880 Sterilit® spray adaptor

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.





GB870R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with Jacobs chuck 0.5 - 7.4 mm,
not for use with acetabulum reamers

0 - 320 rpm; 5 Nm;  3.2 mm;
105 x 29 mm;  230 g

consisting parts:

key GA031R



GB871R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with large AO chuck,
not for use with acetabulum reamers

0 - 320 rpm; 5 Nm;  3.2 mm;
88 x 29 mm;  170 g




GB872R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with Hudson/Zimmer chuck,
not for use with acetabulum reamers

0 - 320 rpm; 5 Nm;  3.2 mm;
86 x 25 mm;  150 g





SMALL DRILL ATTACHMENTS

GB881R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm

0 – 1,250 rpm; 1.5 Nm;  3.2 mm;
89 x 29 mm;  163 g

consisting parts:
key GA031R





GB882R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 4 mm

0 – 1,250 rpm; 1.5 Nm;  3.2 mm;
76 x 22 mm;  100 g

consisting parts:
key GA062R





GB884R

GB716R

DRILL ATTACHMENT

with small A0 chuck

0 – 1,250 rpm; 1.5 Nm;  2.8 mm;
60 x 22 mm;  66 g





GB886R

GB716R

DRILL ATTACHMENT

with AESCULAP® hexagonal chuck

0 – 1,250 rpm; 1.5 Nm;  3.2 mm;
57 x 24 mm;  67 g





GB887R

GB716R

DRILL ATTACHMENT

with Dental shank

0 – 1,250 rpm; 1.5 Nm;  1.8 mm;
56 x 22 mm;  50 g



AESCULAP® Acculan 4

SMALL DRILL ATTACHMENTS

GB878R

GB716R

ADAPTOR

for Synthes radiolucent drive 511.300

0 – 1,250 rpm; 1.5 Nm;
111 x 27 mm; 181 g



GB894R

GB716R

K-WIRE ATTACHMENT

for K-wires from 0.6 – 3.2 mm

0 – 1,250 rpm; 3.2 mm;
92 x 27 x 124 mm; 179 g



GB891R

GB716R

SAGITTAL SAW ATTACHMENT

for "L-sagittal" sawblades with a
maximum usable length of 50 mm

0 – 15,750 rpm; 74 x 33.5 mm;
4.5°; 179 g



GB892R

GB716R

RECIPROCATING SAW ATTACHMENT



for reciprocating saw blades

0 – 15,750 rpm; 135 x 24.5 mm;
166 g; 2.5 mm





GB896R**GB716R****SCREW ATTACHMENT**

with small AO chuck

0 - 165 rpm; 10 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  235 g**GB897R****GB716R****SCREW ATTACHMENT**

with AESCULAP® hexagonal chuck

0 - 165 rpm; 10 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  231 g

AESCULAP® Acculan 4

CHARGER & BATTERIES

GA320

Acculan 4 CHARGER

for Acculan Li-Ion and NiMH batteries GA346, GA676; 4 charging bays; automatic battery diagnosis function; power cord not included

100 – 240 V AC; 50 / 60 Hz;
410 x 120 x 160 mm; 2,300 g



GA346

Acculan 4 LI-ION BATTERY SHORT

only for GA340, GA341 and GA344;
with integrated electronic motor control 9.9 V; 10.89 Wh

120 x 43 x 50 mm; 265 g



GA676

Acculan NIMH BATTERY

only for GA330, GA331, GA334, GA336, GA337; with integrated electronic motor control 9.6 V; 1.95 Ah

140 x 43 x 50 mm; 425 g



GA679

GB489R

Acculan BATTERY REMOVAL DEVICE



GA675

GB487R

Acculan LID

90 g



GA678

GB488R

Acculan STERILE FUNNEL



AESCULAP® Acculan 4 Electro

LOOK. FEEL. EXPERIENCE.



CONSOLE & BATTERIES

GA877

Acculan Electro CONSOLE
power cord not included
100 – 240 V AC; 50 / 60 Hz;
255 x 217 x 90 mm; 2,930 g



GA878 GB486R

Acculan Electro CONNECTION CABLE
for use with Acculan 4 handpieces
3.0 m; 410 g



GA875 GB485R

Acculan Electro LID



GA876

Acculan Electro ADAPTOR
120 x 49 x 42 mm; 165 g



GA678 GB488R

Acculan STERILE FUNNEL



AESCULAP® Acculan 4 ACCESSORIES

ADAPTORS, KEYS & POWER CORD

GB176R GB467R ZIMMER ADAPTOR tool: Zimmer motor: Hudson/Zimmer only in combination with GB630R, GB656R, GB665R, GB669R or GB872R  3.8 mm	GB628R GB690R SYNTHES DHS/DCS ADAPTOR tool: Synthes DHS/DCS motor: AO large  4.2 mm	GB184R GB690R HARRIS ADAPTOR tool: Harris motor: AO large  3.5 mm
XF457R GB690R DIN ADAPTOR tool: DIN motor: AO large  3.8 mm	GA031R KEY for large Jacobs chuck suitable for GB620R, GB621R, GB667R, GB870R, GB881R	GA062R GB690R KEY for small 4 mm Jacobs chuck suitable for GB639R, GB882R
TE780 POWER CORD Europe,  1.5 m		
TE730 POWER CORD Europe,  5.0 m		
TE734 POWER CORD UK,  5.0 m		
TE735 POWER CORD USA, Canada, Japan,  3.5 m		

AESCULAP® PROCESSING

Sterilit® POWER SYSTEMS



AESCULAP® PROCESSING

Sterilit® POWER SYSTEMS

GB600



Sterilit® POWER SYSTEMS

oil spray; 300 ml



GB600850

Sterilit® POWER SYSTEMS SPRAY ADAPTOR

for ELAN 4 air



GB600860

Sterilit® POWER SYSTEMS SPRAY ADAPTOR

for ELAN 4 electro



GB600870

Sterilit® POWER SYSTEMS SPRAY ADAPTOR

for ELAN 4 MIS handpiece shafts



GB600840

Sterilit® POWER SYSTEMS SPRAY ADAPTOR

for Acculan 4



GB600880

Sterilit® POWER SYSTEMS SPRAY ADAPTOR

for Acculan 4 small drill and ELAN 4 electro drill



AESCULAP® MECHANICAL PROCESSING

- The solution for storage and mechanical reprocessing of our AESCULAP® Power Systems.
- For detailed information please refer to the Instruction for Use.



AESCULAP® MECHANICAL PROCESSING

ECCOS® SETS

ELAN 4

GB261R

ELAN 4 electro ECCOS® SET

 106 mm

consisting parts:

GB073R

GB076R


GB083R

GB085R

GB719R

GB072R

JF212R



GB260R

ELAN 4 air ECCOS® SET

 106 mm

consisting parts:

GB053R

GB056R


GB063R

GB065R

GB719R

GB072R

JF212R



GB262R

ELAN 4 electro ECCOS® SET FOR DRILL GA844

 106 mm

consisting parts:

GB072R

GB073R

GB482R


GB715R

GB716R




























JG340

JG342

JG344












ECCOS®
ELAN 4 electro

<p>GB073R</p> <p>ECCOS® HOLDER for motor cable GA804, GA805, GA806 (set consisting of 2 pieces)</p>  56 mm	<p>GB076R</p> <p>ECCOS® HOLDER for perforator driver GA822</p>  56 mm	<p>GB077R</p> <p>ECCOS® HOLDER for intra motor GA824</p>  56 mm
 	 	 
<p>GB082R</p> <p>ECCOS® HOLDER for 1 handpiece</p>  56 mm	<p>GB083R</p> <p>ECCOS® HOLDER for 3 handpieces</p>  56 mm	<p>GB085R</p> <p>ECCOS® HOLDER for craniotome GA849</p>  56 mm
 	 	 
<p>GB086R</p> <p>ECCOS® HOLDER for micro sagittal saw GA836</p>  56 mm	<p>GB087R</p> <p>ECCOS® HOLDER for reciprocating saws GA832, GA837</p>  56 mm	<p>GB088R</p> <p>ECCOS® HOLDER for saws GA831, GA833</p>  56 mm
 	 	 

AESULAP[®] MECHANICAL PROCESSING

ECCOS[®]

ELAN 4 electro

GB482R	GB715R	GB716R
<p>ECCOS[®] HOLDER</p> <p>for ELAN 4 electro drill GA844</p> <p> 106 mm</p>	<p>ECCOS[®] HOLDER</p> <p>for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p> <p> 106 mm</p>	<p>ECCOS[®] HOLDER</p> <p>for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p> <p> 106 mm</p>
 	 	 

ECCOS® & RINSING DEVICES


ELAN 4 electro

GB084R

ECCOS® HOLDER

for MIS handpiece GA860

 56 mm



GB723R

ECCOS® HOLDER

for MIS handpiece shafts GB920R - GB927R


 56 mm



GB692R

RINSING DEVICE


for applied parts, incl. ceramic filter GB693R
(not suitable for autoclaving)



GB679R

RINSING DEVICE

for 4 MIS shafts, incl. ceramic filter GB693R
(not suitable for autoclaving)



GB698R

RINSING ADAPTOR

for applied parts
(for manual cleaning)



GB699R

RINSING ADAPTOR

for MIS shafts
(for manual cleaning)



AESCLAP[®] MECHANICAL PROCESSING

ECCOS[®]

ELAN 4 air

GB053R

ECCOS[®] HOLDER

for motor hose GA705R, GA706R,
GA707R

 56 mm



GB056R

ECCOS[®] HOLDER

for perforator driver GA722

 56 mm




GB057R

ECCOS[®] HOLDER

for intra motor GA724

 56 mm




GB062R

ECCOS[®] HOLDER

for 1 handpiece

 56 mm




GB063R

ECCOS[®] HOLDER

for 3 handpieces

 56 mm



GB065R

ECCOS[®] HOLDER

for craniotome GA749

 56 mm



GB066R

ECCOS[®] HOLDER

for micro sagittal saw GA736

 56 mm



GB067R

ECCOS[®] HOLDER

for reciprocating saws GA732, GA737

 56 mm



GB068R

ECCOS[®] HOLDER

for saws GA731, GA733

 56 mm



ECCOS® & RINSING DEVICES


ELAN 4 air

GB064R

ECCOS® HOLDER

for MIS handpiece GA760

 56 mm



GB723R

ECCOS® HOLDER

for 2 MIS shafts GB920R - GB927R


 56 mm



GB691R

RINSING DEVICE


for applied parts, incl. ceramic filter GB693R
(not suitable for autoclaving)



GB679R

RINSING DEVICE

for 4 MIS shafts, incl. ceramic filter GB693R
(not suitable for autoclaving)



GB697R

RINSING ADAPTOR

for applied parts
(for manual cleaning)



GB699R

RINSING ADAPTOR




for MIS shafts
(for manual cleaning)



AESCULAP® MECHANICAL PROCESSING

ECCOS®

ELAN 4 ACCESSORIES

GB718R ECCOS® HOLDER for 12 1- and 2-ring burrs  56 mm	GB719R ECCOS® HOLDER for 3 craniotome attachments  56 mm	GB720R ECCOS® UNIVERSAL HOLDER for saw blades  56 mm
 	 	 
GB729R ECCOS® HOLDER for 2 reciprocating saw shafts  56 mm	GB730R ECCOS® HOLDER for 2 transversal saw shafts  56 mm	GB072R ECCOS® CABLE/HOSE POS. GUIDE for cable and hoses (set of 4 pieces)  56 mm
 	 	 
GB727R ECCOS® HOLDER for handpiece holding-sleeves GB861P, GB862P, GB863P  56 mm	GB572R ECCOS® FIXATION for GB300R - GB308R  56 mm	
 	 	

ECCOS® SETS

Acculan 4

GB257R

ECCOS® SET

for Acculan 4 small drill GA344

 106 mm

consisting parts:

GB482R

GB487R

GB488R


GB715R

GB716R

JG340

JG342

JG344



GB243R

ECCOS® SET

for 2 Acculan 4 machines

 106 mm

consisting parts:

GB481R

GB487R

GB488R

GB495R

GB496R


GB497R

GB690R

JG340

JG342

JG344



GB256R

ECCOS® SET

for Acculan 4 dermatome

 106 mm


consisting parts:

GB487R

GB488R

GB498R



















JJ464R



AESULAP[®] MECHANICAL PROCESSING

ECCOS[®]

Acculan 4

<p>GB482R</p> <p>ECCOS[®] HOLDER for Acculan 4 small drill GA344</p>  106 mm	<p>GB715R</p> <p>ECCOS[®] HOLDER for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm	<p>GB716R</p> <p>ECCOS[®] HOLDER for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm
 	 	 
<p>GB495R</p> <p>ECCOS[®] HOLDER for GA330, GA331, GA334</p>  106 mm	<p>GB496R</p> <p>ECCOS[®] HOLDER for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R, GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R</p>  106 mm	<p>GB497R</p> <p>ECCOS[®] HOLDER for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R, GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R</p>  106 mm
 	 	 

GB487R**ECCOS® HOLDER**

for GA675

 56 mm
**GB488R****ECCOS® HOLDER**


for GA678

 76 mm
**GB489R****ECCOS® HOLDER**

for GA679

 106 mm
**GB498R****ECCOS® HOLDER**

for GA336, GA337, GA340, GA341

 106 mm
**GB686R****ECCOS® HOLDER**

for GB436R

 56 mm





AESULAP[®] MECHANICAL PROCESSING

ECCOS[®] HOLDER

Acculan 4


GB485R ECCOS[®] HOLDER for GA875  56 mm	GB486R ECCOS[®] HOLDER for GA878  56 mm	GB688R ECCOS[®] HOLDER for BA720R  106 mm
 	 	 
GB481R ECCOS[®] HOLDER for saw blades and other accessories  56 mm	GB690R ECCOS[®] HOLDER 5 mm for drilling wire guide GA672207, GA344244	GB459R ECCOS[®] MOUNTING KIT
 	 	 

ECCOS® HOLDER
Universal Tools


JG303R		JG342 / JG344
REPLACEMENT SCREWS-SET		SILICONE CLAMP
		8.5 mm / 12 mm screw JG303R to be ordered separately
		

AESCULAP® POWER SYSTEMS


TOOL ATTACHMENTS OVERVIEW




AO LARGE




AO SMALL




DENTAL




DIN




HARRIS




HUDSON




HEXAGON



SYNTHES
DHS/DCS














TRINKLE



ZIMMER

SIGNS & SYMBOLS


	<p>LATEX-FREE</p> <p>All products listed in this catalogue are latex-free</p>
	<p>ARTICLE CODE ECCOS® FIXATION</p>
	<p>FORWARD / REVERSE DRIVE</p>
	<p>FORWARD DRIVE</p>
	<p>OSCILLATION FREQUENCY</p>
	<p>OSCILLATON ANGLE</p>
	<p>STROKE</p>
	<p>CANNULATION</p>
	<p>CABLE LENGTH</p>
	<p>WEIGHT</p>
	<p>PACKAGING UNIT</p>
	<p>SINGLE USE</p>
	<p>REQUIRED HEIGHT OF TRAY</p>
	<p>STERILE</p>

AESCULAP® POWER SYSTEMS

BROCHURES

Main Catalogues

084002 Power Systems



017599 Burs & Blades



ELAN 4

071602 ELAN 4



071511 Mechanical Cleaning



012002 Manual Cleaning




Acculan 4


001802 Acculan 4




001902 Acculan 4 Dermotome




002002 Acculan 4 Electro





002102 Acculan 4 Sternum Saws




INDEX

Art. No.	
BA720R	42
GA031R	52
GA062R	52
GA258SU	33
GA259SU	33
GA261SU	33
GA262SU	33
GA263SU	33
GA264SU	33
GA265SU	33
GA266SU	33
GA267SU	33
GA268SU	33
GA269SU	33
GA320	48
GA330	36
GA331	40
GA334	41
GA336	40
GA337	40
GA340	42
GA341	42
GA344	44
GA346	48
GA395SU	33
GA675	48
GA676	48
GA678	48, 51
GA679	48
GA702R	23
GA703R	23
GA705R	24
GA706R	24
GA707R	24
GA708	24
GA710R	23
GA712R	23
GA713R	23
GA715R	23
GA722	24
GA724	31
GA731	31
GA732	31
GA733	31
GA736	31
GA737	31
GA749	25
GA760	29
GA761	26
GA762	26
GA763	26
GA764	26
GA765	27

Art. No.	
GA766	27
GA767	27
GA768	27
GA769	27
GA800	7
GA804	8
GA805	8
GA806	7
GA808	7
GA810	7
GA822	8
GA824	8
GA831	15
GA832	15
GA833	15
GA836	15
GA837	15
GA844	17
GA849	9
GA860	13
GA861	10
GA862	10
GA863	10
GA864	10
GA865	11
GA866	11
GA867	11
GA868	11
GA869	11
GA875	51
GA876	51
GA877	51
GA878	51
GB053R	60
GB056R	60
GB057R	60
GB062R	60
GB063R	60
GB064R	61
GB065R	60
GB066R	60
GB067R	60
GB068R	60
GB072R	62
GB073R	57
GB076R	57
GB077R	57
GB082R	57
GB083R	57
GB084R	59
GB085R	57
GB086R	57
GB087R	57

Art. No.	
GB088R	57
GB176R	52
GB184R	52
GB243R	63
GB256R	63
GB257R	63
GB260R	56
GB261R	56
GB262R	56
GB436R	41
GB459R	66
GB481R	66
GB482R	58, 64
GB485R	66
GB486R	66
GB487R	65
GB488R	65
GB489R	65
GB495R	64
GB496R	64
GB497R	64
GB498R	65
GB572R	62
GB600	54
GB600840	54
GB600850	54
GB600860	54
GB600870	54
GB600880	54
GB620R	36
GB621R	36
GB623R	36
GB628R	52
GB639R	37
GB641R	39
GB645R	37
GB657R	38
GB660R	39
GB663R	37
GB664R	37
GB665R	37
GB667R	38
GB668R	38
GB669R	38
GB670R	38
GB679R	59, 61
GB686R	65
GB688R	66
GB690R	66
GB691R	61
GB692R	59
GB697R	61
GB698R	59

Art. No.	
GB699R	59, 61
GB715R	58, 64
GB716R	58, 64
GB718R	62
GB719R	62
GB720R	62
GB723R	59, 61
GB727R	62
GB729R	62
GB730R	62
GB796SU	33
GB797SU	33
GB861P	32
GB862P	32
GB863P	32
GB866R	32
GB870R	17, 44
GB871R	17, 44
GB872R	17, 44
GB878R	19, 46
GB881R	18, 45
GB882R	18, 45
GB884R	18, 45
GB886R	18, 45
GB887R	18, 45
GB891R	19, 46
GB892R	19, 46
GB894R	19, 46
GB896R	20, 47
GB897R	20, 47
GB920R	13, 29
GB921R	13, 29
GB922R	13, 29
GB925R	13, 29
GB926R	13, 29
GB927R	13, 29
GB941R	9, 25
GB942R	9, 25
GB943R	9, 25
GB945R	9, 25
GB947R	9, 25
JG303R	67
JG342	67
JG344	67
TE730	32, 52
TE734	32, 52
TE735	32, 52
TE780	32, 52
XF457R	52

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "AESCULAP" and the product trademarks "ECCOS" and "Sterilit" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.