

Bone Wax

Bone Wax Stick



Biosurgicals

Bone Wax

Bone Wax Stick



Classic Bar or New Stick

New, easier, stick format ...

- The new **Bone Wax Stick** format allows spreading of the product onto the bleeding area in a comfortable and clean way.

... or Classic bars

- The classic, widely used **Bone Wax** bars
- Can be moulded and kneaded only **hand-warm** and applied to the bone using a spatula.

Ordering Information

Description	Code
12 Stick set	102 9755
24 Bars set	102 9754

Features


Sterile mixture of natural beeswax and vaseline

- ▶ Presentation in a new **application stick** or the classic **kneading bars**
- ▶ A small amount is sufficient to reach haemostasis
- ▶ Easy mechanical stilling of bleeding at bones
- ▶ Total hemostasis by tamponade effect
- ▶ Achieve mechanical occlusion of intraosseous vessels

Indicated for mechanical stilling of bleeding at bones in:

Indications

- ▶ Thoracic surgery (sternum & ribs)
- ▶ Dental, oral and jaw surgery
- ▶ Orthopaedics and traumatology
- ▶ Neurosurgery (trepanation)





Distributed by

B. Braun Surgical SA | Carretera de Terrassa, 121 | 08191 Rubí | Spain
Phone +34 93 5 86 62 00 | Fax +34 93 6 99 73 03 | www.bbraun.es

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 07461 95-0 | Fax +49 07461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products.

No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.