

Bone Wax

For mechanical haemostasis in bone injuries

Description

Bone Wax is a sterile mixture of beeswax (70%) and vaseline (30%), which can be moulded and applicated in a "hand-warm" condition.


Indication

Bone Wax is used for the mechanical haemostasis in bones:

- In neurosurgery
- In orthopaedics and traumatology
- In thoracic surgery
- In dental, oral and jaw surgery

Advantages

- Soft, easy to form and to apply
- Easy to handle due to practical, transparent double peel-packs



Content / Box

24 pieces each 2.5 gr.

Cat. N°

102 9754

Lyostypt®

Gelita Tampon

Bone Wax