




NEUROSURGERY

AESCULAP[®] SENSATION MICRO INSTRUMENTS

OUR TIP FOR THE TOP

AESCULAP® SENSATION MICRO INSTRUMENTS

OUR TIP FOR THE TOP



This pictogram shows the type of instrument.

Colored

Color-coded profile handles enable a clear identification of the type of instrument.



Noir® scissors



Scissors






Tumor grasping
forceps




Forceps



Needle holders



Angled bayonet shape
For enhanced sight lines.
It moves the surgeons hand further out of the field of vision while working under the microscope.



AESCULAP® SENSATION MICRO INSTRUMENTS

Noir® – OUR TIP FOR THE TOP

The well known AESCULAP® Noir®-coating offers the advantage that disturbing light reflections can be reduced.



Noir® Scissors,
upwards curved




Serrated blades reduce
the risk of tissue slipping
out of the jaws



RED color-coded
profile handle
for Noir® scissors





SCISSORS



Jaws


sharp/sharp 


	FM146B	FM147B	FM148B	FM149B
Working length	60 mm 2 1/8"	70 mm 2 3/4"	90 mm 3 1/2"	120 mm 4 3/4"
Total length	185 mm 7 1/3"	195 mm 7 3/4"	215 mm 8 1/2"	245 mm 9 3/4"

1/2

AESCULAP® SENSATION MICRO INSTRUMENTS

SILVER – OUR TIP FOR THE TOP









SCISSORS

Scissors,
downwards curved

Scissors,
angled



1/2

Jaws

sharp/sharp	
blunt/blunt	


	FM163R			
	FM164R	FM167R	FM168R	FM169R
Working length	90 mm 3 1/2"	90 mm 3 1/2"	90 mm 3 1/2"	90 mm 3 1/2"
Total length	220 mm 8 1/2"	220 mm 8 1/2"	220 mm 8 1/2"	220 mm 8 1/2"

AESCULAP® SENSATION MICRO INSTRUMENTS

SILVER – OUR TIP FOR THE TOP



Scissors,
straight





Jaws

sharp/sharp	
sharp/blunt	
blunt/blunt	

	FM120R	FM122R	FM124R	FM159R
	FM130R	FM132R	FM134R	-
	FM140R	FM142R	FM144R	FM160R
Working length	70 mm 2 3/4"	90 mm 3 1/2"	120 mm 4 3/4"	135 mm 5 1/3"
Total length	195 mm 7 3/4"	220 mm 8 1/2"	245 mm 9 3/4"	265 mm 10 1/2"



Scissors,
upwards curved



Jaws


sharp/sharp	
sharp/blunt	
blunt/blunt	

	FM121R	FM123R	FM125R	FM161R
	FM131R	FM133R	FM135R	-
	FM141R	FM143R	FM145R	FM162R
Working length	70 mm 2 3/4"	90 mm 3 1/2"	120 mm 4 3/4"	135 mm 5 1/3"
Total length	195 mm 7 3/4"	220 mm 8 1/2"	245 mm 9 3/4"	265 mm 10 1/2"


1/2

AESCULAP® SENSATION MICRO INSTRUMENTS

SILVER – OUR TIP FOR THE TOP



GRAY color-coded profile
handle for grasping forceps




Forceps with teeth for grasping
and holding of tissue and for
soft lifting of fine structures



Serrated ring tips for tumor
removal



TISSUE AND TUMOR GRASPING FORCEPS



Jaws			
 1 x 2 teeth	Tissue forceps	FM174R	
 2.5 mm	Tumor grasping forceps	-	FM176R FM178R
 3.5 mm	Tumor grasping forceps	-	FM177R FM179R
Working length		90 mm 3 1/2"	90 mm 3 1/2" 120 mm 4 3/4"
Total length		210 mm 8"	210 mm 8" 245 mm 9 3/4"

AESCULAP® SENSATION MICRO INSTRUMENTS

SILVER – OUR TIP FOR THE TOP



BLACK color-coded profile
handle for forceps



Grasping of fine structures



Pin to avoid any scissoring



FINE FORCEPS



1/2

Jaws

0.5 mm

0.9 mm


Working length

Total length

FM150R	FM153R	FM156R
FM151R	FM154R	FM157R
70 mm 2 3/4"	90 mm 3 1/2"	120 mm 4 3/4"
195 mm 7 3/4"	215 mm 8 1/2"	245 mm 9 3/4"

AESCULAP® SENSATION MICRO INSTRUMENTS


SILVER – OUR TIP FOR THE TOP



GREEN color-coded profile
handle for needle holders



Curved and smooth needle
holder tip for needle sizes 7/0
and smaller



For storing SENSATION micro instruments we recommend FM665R perforated basket. Up to 14 instruments can be stored.




JK701P

3/4 container (basic version) suitable for one FM665R (without illustration).

FM665R

Perforated basket with instrument racks and lid for storage of SENSATION micro instruments (without instruments).



Jaws			
0.5 mm			
			
	FM170R	FM173R	FM175R
Working length	70 mm 2 3/4"	90 mm 3 1/2"	120 mm 4 3/4"
Total length	198 mm 7 3/4"	217 mm 8 1/2"	247 mm 9 3/4"

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "AESCULAP" and the product trademark "Noir" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.