


# B. Braun Vascular Systems

## Uni-Graft® K DV Patch

Impregnated Polyester Patch  
for Reconstructive Interventions






# Uni-Graft® K DV Patch

## Zero porosity patch material for arterial reconstructions

B. Braun has developed an impregnated vascular patch from proven polyester. The Uni-Graft® K DV Patch meets the requirements of a synthetic patch for vascular surgery:

- good handling properties
- high tensile strength and elasticity
- good incorporation and compatibility
- long-term stability
- low thrombogenicity

Uni-Graft® K DV Patch is a warp-knitted double-velour patch with an impregnation of absorbable modified gelatin. Uni-Graft® K DV Patch is absolutely blood-tight at the time of implantation.



The unique modified gelatin impregnation reduces suture hole bleeding.



Dimension	Content	Article Number
1 x 7 cm	1 piece	110 7401
	2 pieces	110 7404
	5 pieces	110 7407
2 x 9 cm	1 piece	110 7400
	2 pieces	110 7403
	5 pieces	110 7406
3 x 4 cm	1 piece	110 7402
	2 pieces	110 7405
	5 pieces	110 7408

Uni-Graft® K DV Patch combines a number of significant advantages:

- minimised suture hole bleedings
- no tearing at the suture line due to the high tensile strength
- good moulding properties and flexibility
- excellent tissue incorporation

For Product information please contact:

**B. Braun Melsungen AG**

Vascular Systems

Sieversufer 8

D-1235 Berlin

Phone +49 30 689897-0

Fax +49 30 689897-30

**B | BRAUN**  
SHARING EXPERTISE

**AESCLAP AG & CO. KG**

Am Aesculap-Platz  
78532 Tuttlingen

Phone +49 7461 95-0

Fax +49 7461 95-2600

[www.aesculap.de](http://www.aesculap.de)

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.