

B. Braun Direct Dispense Packaging

A functional and modern packaging design



B. Braun Direct Dispense Packaging

Clearly labelled information so that the suture material can be identified rapidly and reliably.


The most important product data is printed twice on the pack (horizontally and vertically).



Smooth and transparent outer peel pack.



Tear off aluminium opening system.




Reduced waste
and consequently
secure disposal.
Only one piece is
involved in the
disposal of the
packaging
material.



A very rapid and
reliable, one at a
time removal pack.
Pre-cuts can be
removed
individually. No
tangling, no
knotting, right
down to the last
suture.



Excellent needle
positioning for a
quick and direct
gripping with the
needle holder.
A very smooth and
easy removal is
guaranteed.



Tear along arrow
marked, with a
single movement.
Pull the tab
downwards.
No force is needed
for the opening.



Box Label

Suture name: Safil® Quick

Thread description: Synthetic, fast absorbable surgical suture, polyglycolic acid, coated, braided, UNDYED

USP size: 2/0

Metric size: metric 3

Thread length cm/inches: 70 cm (28")

Responsible manufacturer: B|BRAUN AEscuLAP
AESCuLAP AG & CO. KG
D-78532 Tuttlingen

Sterilization method: Sterile R: Gamma irradiation
Sterile EO: Ethylene oxide

"Do not reuse"

Contains instructions for use

Recyclable

CE Mark Registration: 0123

Product code: C0046042

Units per box: 36 units
G=12 units
B=24 units
C=36 units

Batch number: LOT 1/0253

Expiry date: 2005-06

Needle description: HR 26
1/2 c 26 mm
Taper

Easyslide

300 series stainless steel needle

Absorbables: Safil®, Safil® Quick, Monosyn®, Softcat® Plain, Softcat® Chrom

Non-Absorbables: Premilene® Securex, Dafilon® Nylon, Supramid, Synthofil® PremiCron®, Dagrofil®, Miralene®, Silekam® Virgin Silk, Steelex Sternum Set, Steelex Electrode Set, Linatrix®

B|BRAUN
AESCuLAP

AESCuLAP AG & CO KG
Am Aesculap-Platz
D-78532 Tuttlingen
Phone +49 (7461) 95-0
Fax +49 (7461) 95-2600
Internet <http://www.aesculap.de>